

# INNILOCK MONOAXIAL SCREW 7.5MM ( FOR 6MM ROD )



CAT. NO.	METAL	INNILOCK Monoaxial Screw 7.5mm ( For 6mm rod )
643.730	TT	INNILOCK Monoaxial Screw - Ø 7.5 mm x 30 mm
643.735	TT	INNILOCK Monoaxial Screw - Ø 7.5 mm x 35 mm
643.740	TT	INNILOCK Monoaxial Screw - Ø 7.5 mm x 40 mm
643.745	TT	INNILOCK Monoaxial Screw - Ø 7.5 mm x 45 mm
643.750	TT	INNILOCK Monoaxial Screw - Ø 7.5 mm x 50 mm

## Specification

- Available length 30,35,40,45 & 50mm with 7.5mm dia.

## Indication for use

- The 7.5 MM Monoaxial screw is indicated for degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies).
- Pedicle screw fixation is limited to skeletally mature patients.
- The 7.5 MM Monoaxial Screw is a non-cervical spinal fixation device intended for use as a posterior pedicle screw.