

LATERAL MASS / POSTERIOR CERVICAL POLY SCREWS 4.0 MM



SR. NO.	CAT. NO. TT	METAL	LATERAL MASS / POSTERIOR CERVICAL POLY SCREWS 4.0 mm
1	695.110	TT	Lateral Mass / Posterior cervical Poly screws 4.0 mm / 10 mm
2	695.112	TT	Lateral Mass / Posterior cervical Poly screws 4.0 mm / 12 mm
3	695.114	TT	Lateral Mass / Posterior cervical Poly screws 4.0 mm / 14 mm
4	695.116	TT	Lateral Mass / Posterior cervical Poly screws 4.0 mm / 16 mm
5	695.118	TT	Lateral Mass / Posterior cervical Poly screws 4.0 mm / 18 mm
6	695.120	TT	Lateral Mass / Posterior cervical Poly screws 4.0 mm / 20 mm
7	695.122	TT	Lateral Mass / Posterior cervical Poly screws 4.0 mm / 22 mm
8	695.124	TT	Lateral Mass / Posterior cervical Poly screws 4.0 mm / 24 mm
9	695.126	TT	Lateral Mass / Posterior cervical Poly screws 4.0 mm / 26 mm
10	695.128	TT	Lateral Mass / Posterior cervical Poly screws 4.0 mm / 28 mm
11	695.130	TT	Lateral Mass / Posterior cervical Poly screws 4.0 mm / 30 mm
12	695.136	TT	Lateral Mass / Posterior cervical Poly screws 4.0 mm / 36 mm

Specification

- Lateral Mass / Posterior cervical Poly screws Available in Titanium.
- Available length from 10,12,14,16,18,20,22,24,26,28,30 & 36mm with 4.0mm diameter.

Indication for use

- Lateral Mass Poly axial screw is used for stabilization of the occipitocervical junction.

