

INNILOCK MONOAXIAL SCREW 4.5MM (FOR 6MM ROD)



CAT. NO.	METAL	INNILOCK Monoaxial Screw 4.5mm (For 6mm rod)
643.025	TT	INNILOCK Monoaxial Screw Ø 4.5 mm x 25 mm
643.030	TT	INNILOCK Monoaxial Screw Ø 4.5 mm x 30 mm
643.035	TT	INNILOCK Monoaxial Screw Ø 4.5 mm x 35 mm
643.040	TT	INNILOCK Monoaxial Screw Ø 4.5 mm x 40 mm
643.045	TT	INNILOCK Monoaxial Screw Ø 4.5 mm x 45 mm

Specification

- Available length 25, 30, 35, 40 & 45mm with 4.5mm dia.

Indication for use

- The 4.5 MM Monoaxial screw is indicated for degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies).
- Pedicle screw fixation is limited to skeletally mature patients.
- The 4.5 MM Monoaxial Screw is a non-cervical spinal fixation device intended for use as a posterior pedicle screw.