

INNILOCK POLYAXIAL SCREW- Ø 5.5 MM (6MM ROD SYSTEM)



SR.NO.	CAT. NO. TT	METAL	INNILOCK POLYAXIAL SCREW - Ø 5.5 mm (6mm rod system)
1	644.13 0	TT	INNILOCK Polyaxial Screw - Ø 5.5 mm x 30 mm
2	644.135	TT	INNILOCK Polyaxial Screw - Ø 5.5 mm x 35mm
3	644.14 0	TT	INNILOCK Polyaxial Screw - Ø 5.5 mm x 40 mm
4	644.145	TT	INNILOCK Polyaxial Screw - Ø 5.5 mm x 45 mm
5	644.15 0	TT	INNILOCK Polyaxial Screw - Ø 5.5 mm x 50 mm

Specification

- It is a Spine Fixation System Intended for use as a pedicle Screw fixation or as an Anterolateral lateral spine fixation system.
- Available length 30,35,40,45 & 50mm.

Indication for use

- The 5.5 MM Polyaxial screw is indicated for degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies). Levels of fixation are for the Thoracic, Lumbar and Sacral spine.