

CANNULATED DRILL TAP WITH T-HANDLE- CALIBRATED



| CAT. NO. | METAL | Cannulated Drill Tap with T-Handle- Calibrated |
|----------|-------|---|
| 429.435 | SS | Cannulated Drill Tap with T-Handle Ø 3.5 mm |
| 429.440 | SS | Cannulated Drill Tap with T-Handle Ø 4.0 mm |
| 429.445 | SS | Cannulated Drill Tap with T-Handle Ø 4.5 mm |
| 429.669 | SS | Cannulated Drill Tap with T-Handle, Calibrated Ø 6.5 mm |
| 429.670 | SS | Cannulated Drill Tap with T-Handle, Calibrated Ø 7.0 mm |

Specification

- Cannulated Drill Tap with T-Handle is available in Stainless Steel.
- Available with different diameter 3.5, 4.0, 4.5, 6.5 & 7mm.

Indication for use

- Cannulated Drill Tap with T-Handle is used in orthopedic surgery for holding taps, drills, reamers and other small tools to be turned by hand.