



CAT. NO.	METAL	DRILL SLEEVE, TAP SLEEVE
1596.005	SS	Double Drill Sleeve Ø 3.5/2.5 mm
1596.006	SS	Double Drill Sleeve Ø 4.0/2.5 mm
1596.007	SS	Double Drill Sleeve Ø 4.5/3.2 mm
1596.008	SS	Double Drill Sleeve Ø 6.5/3.2 mm
1596.009	SS	Double Drill Sleeve Ø 6.5/4.5 mm
1596.065	SS	TAP SLEEVE 3.5MM FOR TAP 3.5MM
1596.070	SS	DRILL SLEEVE Ø 0.8MM FOR 1.2MM CORTICAL SCREW
312.470	SS	Drill Sleeve Pointed (€Pointed Drill Guide)
312.480	SS	Insert Drill Sleeve 4.5/ 3.2mm

## Specification

- It is available in Stainless Steel.

## Indication for use

- Drill guides and sleeves are used with drills and taps to place accurate holes and protect surrounding tissue.

Tap sleeve can also be used as a drill bit guide.

