



CAT. NO.	METAL	Drill Tap with T-Handle
429.015	SS	Drill Tap with T-Handle Ø 1.5 mm
429.022	SS	Drill Tap with T-Handle Ø 2.0 mm
429.024	SS	Drill Tap with T-Handle Ø 2.4 mm
429.027	SS	Drill Tap with T-Handle Ø 2.7 mm
429.035	SS	Drill Tap with T-Handle Ø 3.5 mm (20 TPI)
429.040	SS	Drill Tap with T-Handle Ø 4.0 mm (14 TPI, cancellous)
429.045	SS	Drill Tap with T-Handle Ø 4.5 mm
429.065	SS	Drill Tap with T-Handle Ø 6.5 mm (9 TPI, cancellous)
429.135	SS	Drill Tap with T-Handle Ø 3.5 mm (14 TPI, cancellous)
429.140	SS	Drill Tap with T-Handle Ø 10.5 mm FOR PFNA II NECK SCREWS

Specification

- Drill Tap with T-Handle is available in Stainless Steel.
- Available with different diameter 1.5,2.0,2.4,2.7,3.5,4.0,4.5,6.5 & 10.5mm.

Indication for use

- Drill Tap with T-handle tap is used in orthopaedic surgery for holding taps, drills, reamers and other small tools to be turned by hand.

