



CAT. NO.	METAL	FORCEPS
420.035	SS	Screw Holding Forcep- Ø 3.5 mm
420.045	SS	Screw Holding Forcep- Ø 4.5 mm
420.080	SS	Forcep for Screw Removal 3.5MM
420.090	SS	Forcep for Screw Removal 4.5MM

Specification

- Forceps is available in Stainless steel.
- It is available in 3.5 - 4.5mm dia.

Indication for use

- Forceps is used in orthopedic surgery for holding the screws perfectly.

As well, it is used for picking up bone screws from the rack of screws.

It will also facilitate the user and reduce the surgery time.