



CAT. NO.	METAL	Flexible Reamer Shaft
441.059	SS	FLEXIBLE SHAFT FOR DRILL BIT
443.080	SS	Flexible Reamer Shaft- Ø 8.0 mm
443.090	SS	Flexible Reamer Shaft- Ø 9.0 mm
443.100	SS	Flexible Reamer Shaft- Ø 10.0 mm
443.101	SS	Flexible Reamer Shaft 24"

Specification

- Flexible Reamer Shaft is available in Stainless steel.
- Available Dia 8,9 & 10mm.

Indication for use

- The flexible reamer shaft is intended to clear the medullary canal of bone marrow and debris and to effectively enlarge the medullary canal of long bones for the acceptance of an intramedullary nail.

