



CAT. NO.	METAL	GAUGES
1598.010	SS	DEPTH Gauge FOR 1.2MM SCREWS
1598.014	SS	DEPTH Gauge FOR 1.5MM SCREWS
421.030	SS	Depth Gauge for Ø 2.4/2.7 mm
421.034	SS	Depth Gauge for Ø 3.4 mm
421.035	SS	Depth Gauge for Ø 3.5 mm
421.045	SS	Depth Gauge for Ø 4.5 mm
903.150	SS	Depth Gauge For Interlocking- Ø 4.9mm

Specification

- Gauge is available in Stainless steel.
- It is available in 2.4/2.7, 3.4, 3.5, 4.5 & interlocking 4.9mm dia.

Indication for use

- A Depth Gauge is a specifically designed surgical instrument that is very precise and efficient in measuring depths during orthopedic surgeries. This instrument is primarily used to gauge the depth of the holes drilled in bones during implants and other similar surgical procedures.

