



CAT. NO.	METAL	HAMMER
1856.012	SS	Bone Hammer with Fibre Handle 200 gms.
1856.015	SS	Bone Hammer with Fibre Handle 350 gms.
1856.025	SS	Bone Hammer with Fibre Handle 500 gms.
1856.027	SS	Bone Hammer with Fibre Handle 700 gms.
1856.035	SS	Bone Hammer with S.S. Handle 350gms.
1856.051	SS	Bone Hammer Nylon Faced
904.130	SS	SLOTTED HAMMER FOR HUMERUS NAILS

Specification

- Bone Hammer is available in Stainless steel.
- It is available in 200,350,500 & 700 gms with Fibre Handle.

Indication for use

- Bone Hammer is an orthopedic surgical instrument commonly used by orthopedic surgeons to impact or exert force on osteotomes, chisels, gouges, tampers, and other specially designed instruments.

Bone Hammer has a Fibre handle to provide convenient handgrip.