



CAT. NO.	METAL	LC-LCP Drill Sleeve
1596.104	SS	LC-LCP Drill Sleeve/ Guide Ø 1.0 mm for lcp 1.5mm screws
415.355	SS	LC-LCP Trochantric Sleeve 3.2mm
910.102	SS	LC-LCP Drill Sleeve/ Guide Ø 2.8 mm
910.104	SS	LC-LCP Combined Drill Sleeve Ø 2.5mm / 2.8 mm
910.105	SS	LC-LCP Combined Drill Sleeve Ø 2.0mm / 2.7 mm
910.106	SS	LC-LCP Combined Drill Sleeve Ø 1.8mm / 2.4 mm
910.107	SS	LC-LCP Drill Sleeve/ Guide Ø 1.8 mm
910.108	SS	LC-LCP Drill Sleeve/ Guide Ø 2.0 mm
910.109	SS	DISPOSABLE SLEEVE 2.8MM
910.110	SS	DISPOSABLE SLEEVE 4.3MM
910.111	SS	DISPOSABLE SLEEVE 1.8MM
910.112	SS	DISPOSABLE SLEEVE 2.0MM
910.113	SS	LC-LCP Drill Sleeve/ Guide Ø 1.5 mm for lcp 2.0mm screws
910.128	SS	LC-LCP Drill Sleeve/ Guide Ø 4.3 mm
910.129	SS	LC-LCP Combined Drill Sleeve Ø 3.2mm / 4.3mm
910.131	SS	LCP VAIRABLE ANGLE DRILL SLEEVE 1.8MM
910.132	SS	LCP VAIRABLE ANGLE DRILL SLEEVE 2.8MM

## Specification

- LC-LCP Drill Sleeve is available in Stainless steel.
- It is available with 1.8,2.0,2.5 & 2.8mm Dia.

## Indication for use

- The use of the LC-LCP drill sleeve is mandatory in order to ensure that the locking screw is drilled in the proper perpendicular angle and correctly locked in the plate.

