



CAT. NO.	METAL	Osteotome with Fibre Handle
411.050	SS	Osteotome with Fibre Handle-5 mm, Straight
411.051	SS	Osteotome with Fibre Handle-5 mm, Curved
411.060	SS	Osteotome with Fibre Handle-6 mm, Straight
411.061	SS	Osteotome with Fibre Handle-6 mm, Curved
411.100	SS	Osteotome with Fibre Handle- 10 mm, Straight
411.101	SS	Osteotome with Fibre Handle- 10 mm, Curved
411.120	SS	Osteotome with Fibre Handle- 12 mm, Straight
411.121	SS	Osteotome with Fibre Handle- 12 mm, Curved
411.150	SS	Osteotome with Fibre Handle- 15 mm, Straight
411.151	SS	Osteotome with Fibre Handle- 15 mm, Curved
411.180	SS	Osteotome with Fibre Handle- 18 mm, Straight
411.181	SS	Osteotome with Fibre Handle- 18 mm, Curved
411.200	SS	Osteotome with Fibre Handle- 20 mm, Straight
411.201	SS	Osteotome with Fibre Handle- 20 mm, Curved
411.240	SS	Osteotome with Fibre Handle- 24 mm, Straight
411.250	SS	Osteotome with Fibre Handle- 25 mm, Straight
411.251	SS	Osteotome with Fibre Handle- 25 mm, Curved
411.300	SS	Osteotome with Fibre Handle- 30 mm, Straight
411.301	SS	Osteotome with Fibre Handle- 30 mm, Curved

## Specification

- Osteotome with Fibre Handle is available in Stainless steel.
- It is available in 5,6,10,12,15,18,20,24,25 & 30mm length..

## Indication for use

- The GPC osteotomes with fibre handle is a osteotome with both ends which are sharp. It is used for osteotomy of the bone eg: HTO, corrective osteotomy, modified french osteotomy.