



CAT. NO.	METAL	Osteotome with S.S Handle
1761.006	SS	Osteotome with S.S Handle-6 mm, Straight
1761.008	SS	Osteotome with S.S Handle-8 mm, Straight
1761.012	SS	Osteotome with S.S Handle- 12 mm, Straight
1761.013	SS	Osteotome with S.S Handle- 13 mm, Straight
1761.015	SS	Osteotome with S.S Handle- 15 mm, STRAIGHT
1761.018	SS	Osteotome with S.S Handle- 18 mm, Straight
1761.020	SS	Osteotome with S.S Handle- 20 mm, STRAIGHT
1761.024	SS	Osteotome with S.S Handle- 24 mm, Straight
1761.025	SS	Osteotome with S.S Handle- 25 mm, Straight
412.061	SS	Osteotome with S.S Handle-6 mm, Curved
412.081	SS	Osteotome with S.S Handle-8 mm, Curved
412.121	SS	Osteotome with S.S Handle- 12 mm, Curved
412.151	SS	Osteotome with S.S Handle- 15 mm, Curved
412.181	SS	Osteotome with S.S Handle- 18 mm, Curved
412.201	SS	Osteotome with S.S Handle- 20 mm, Curved
412.241	SS	Osteotome with S.S Handle- 24 mm, Curved
412.300	SS	Osteotome with S.S Handle- 30 mm, Straight
412.301	SS	Osteotome with S.S Handle- 30 mm, Curved

Specification

- Osteotome with S.S Handle is available in Stainless steel.
- It is available in 6,8,12,13,15,18,20,24,25 & 30mm length.

Indication for use

- An osteotome is an instrument used for cutting or preparing bone. osteotomes are similar to a chisel but beveled on both sides. they are used orthopedic surgery.

