



CAT. NO.	METAL	Patella Holding Forcep
1723.617	SS	Patella Holding Forcep- 175mm
406.150	SS	Patella Holding Forcep- 150MM
406.190	SS	Patella Holding Forcep- 190mm
406.210	SS	Patella Holding Forcep- 210mm
406.240	SS	Patella Holding Forcep- 240mm

## Specification

- It is available in 175,150,190,210 & 240mm length.

## Indication for use

- Patella forceps are ideally designed for use in knee surgeries. procedures such as patellar realignment or treatments for patellar fractures and dislocations can make use of this instrument.

