



CAT. NO.	METAL	REAMER
433.006	SS	K. Nail Reamer Ø 6 mm
433.007	SS	K. Nail Reamer Ø 7 mm (20")
433.008	SS	K. Nail Reamer Ø 8 mm
433.009	SS	K. Nail Reamer Ø 9 mm
433.010	SS	K. Nail Reamer Ø 10 mm
433.011	SS	K. Nail Reamer Ø 11 mm
433.012	SS	K. Nail Reamer Ø 12 mm
433.013	SS	K. Nail Reamer Ø 13 mm
433.014	SS	K. Nail Reamer Ø 14 mm
433.107	SS	K. Nail Reamer Ø 7 mm (22")
433.200	SS	Square Nail Reamer Ø 2.0 mm
433.225	SS	Square Nail Reamer Ø 2.5 mm
433.230	SS	Square Nail Reamer Ø 3.0 mm
433.235	SS	Square Nail Reamer Ø 3.5 mm
433.240	SS	Square Nail Reamer Ø 4.0 mm
909.155	SS	Cannulated Hip pin Reamer 6.2 mm
909.160	SS	Cannulated Femoral Neck Screw Reamer 8mm

Specification

- Reamer is available in Stainless steel.
- K.Nail Reamer available with 6-14mm dia & square nail available with 2-4mm dia.

Indication for use

- The reamer system is intended to clear the medullary canal of bone marrow and debris and to effectively enlarge the medullary canal of long bones for the acceptance of an intramedullary nail.