



| CAT.NO. | METAL | RETRACTORS |
|---------|-------|--|
| 435.005 | SS | SelfRetaining Retractor- Small |
| 435.010 | SS | Self Retaining Retractor- Medium |
| 435.015 | SS | Self Retaining Retractor- Large |
| 435.025 | SS | SelfRetaining Retractor- Ex. Large |
| 435.050 | SS | Nerve Root Retractor (Set of 3) |
| 435.060 | SS | Nerve Root Retractor- Straight |
| 435.065 | SS | Nerve Root Retractor- 45 deg |
| 435.070 | SS | Nerve Root Retractor- 90 deg |
| 435.080 | SS | Bristow Retractor |
| 435.090 | SS | Dura Retractor |
| 435.100 | SS | Charnley Retractor with Weight & Chain |
| 435.200 | SS | Gelpi Retractor |
| 659.104 | SS | Bilateral Spine Retractor |

Specification

- Retractors are available in Stainless steel.
- Self Retaining Retractors are available in small, medium, large & extra large size.

Indication for use

- RETRACTORS is instrument for medical. A retractor is a surgical instrument by which a surgeon can either actively separate the edges of a surgical incision or wound, or can hold back underlying organs and tissues, so that body parts under the incision may be accessed.