



CAT. NO.	METAL	Rongeur
438.020	SS	Up Cutting Rongeur- Ø 2.0 mm
438.025	SS	Up Cutting Rongeur- Ø 2.5 mm
438.030	SS	Up Cutting Rongeur- Ø 3.0 mm
438.040	SS	Up Cutting Rongeur- Ø 4.0 mm
438.120	SS	Down Cutting Rongeur- Ø 2.0 mm
438.125	SS	Down Cutting Rongeur- Ø 2.5 mm
438.130	SS	Down Cutting Rongeur- Ø 3.0 mm
438.140	SS	Down Cutting Rongeur- Ø 4.0 mm

Specification

- Rongeur is available in Stainless steel.
- Available Diameter 2.0,2.5,3.0 & 4.0mm.

Indication for use

- A Rongeur can be used to open a window in bone, often in the skull, in order to access tissue underneath. They are used in neurosurgery, podiatric surgery, maxillofacial surgery, and orthopedic surgery to expose areas for operation.