



CAT. NO.	METAL	Self Centering Bone Holding Forcecp
405.150	SS	Self Centering Bone Holding Forcecp- 150 mm
405.170	SS	Self Centering Bone Holding Forcecp- 170 mm
405.175	SS	Self Centering Bone Holding Forcecp- 175 mm
405.190	SS	Self Centering Bone Holding Forcecp- 190 mm
405.200	SS	Self Centering Bone Holding Forcecp- 200 mm
405.240	SS	Self Centering Bone Holding Forcecp- 240 mm
405.260	SS	Self Centering Bone Holding Forcecp- 260 mm
405.270	SS	Self Centering Bone Holding Forcecp- 270 mm
405.280	SS	Self Centering Bone Holding Forcecp- 280 mm

Specification

- Self Centering Bone Holding Forcecp is available in Stainless steel.
- It is available in 150-280mm length.

Indication for use

- Self Centering Bone Holding Forceps is an orthopedic surgical tool that surgeons commonly use to clamp, hold and move fractured bone fragments, in order to reposition them into their normal alignment.

