



CAT. NO.	ADAPTION PLATE 20 HOLES
706.600	ADAPTION PLATE 20 HOLES THICKNESS 0.6 MM
707.080	ADAPTION PLATE 20 HOLES THICKNESS 0.7 MM
707.650	ADAPTION PLATE 20 HOLES THICKNESS 0.8 MM

Specification

- Plate is available in Titanium.

Indication for use

- An Adaption Plate with 20 holes is used in oral and maxillofacial surgery. It is primarily employed for stabilizing fractures, reconstructing bone, or supporting grafts, especially in cases where a longer span or multiple points of fixation are required.