



CAT. NO.	DOUBLE -Y-PLATE 6 HOLES
706.615	DOUBLE -Y-PLATE 6 HOLES THICKNESS 0.6 MM
707.096	DOUBLE -Y-PLATE 6 HOLES THICKNESS 0.7 MM

Specification

- Plate is available in Titanium.

Indication for use

- A Double-Y Plate with 6 holes is a specialized fixation plate used in oral and maxillofacial surgery, particularly designed for stabilizing mandibular condylar fractures or condylar neck fractures. Its unique "double Y" shape provides multidirectional fixation and better distribution of forces across the fracture site.