



CAT. NO.	FRAME PLATE 4 HOLES (05X10 MM)
706.620	FRAME PLATE 4 HOLES (05X10 MM) THICKNESS 0.6 MM
707.101	FRAME PLATE 4 HOLES (05X10 MM) THICKNESS 0.7 MM

Specification

- Plate is available in Titanium.

Indication for use

- A Frame Plate with 4 holes (0.5 x 10 mm) is a small, low-profile fixation plate used in oral and maxillofacial surgery, especially for the stabilization of small bone segments such as in alveolar ridge fractures, paediatric cases, or delicate areas like the orbital rim or nasal bones.