



CAT. NO.	FRAME PLATE 4 HOLES (10X10 MM)
706.625	FRAME PLATE 4 HOLES (10X10 MM) THICKNESS 0.6 MM
707.106	FRAME PLATE 4 HOLES (10X10 MM) THICKNESS 0.7 MM

Specification

- Plate is available in Titanium.

Indication for use

- A Frame Plate with 4 holes is a small, low-profile fixation plate used in oral and maxillofacial surgery, especially for the stabilization of small bone segments such as in alveolar ridge fractures, paediatric cases, or delicate areas like the orbital rim or nasal bones.

