



CAT NO.	H-PLATE 11 HOLES
706.501	H-PLATE 11 HOLES THICKNESS 0.6MM
707.001	H-PLATE 11 HOLES THICKNESS 0.7MM
707.601	H-PLATE 11 HOLES THICKNESS 0.8MM

Specification

- Plate is available in Titanium.

Indication for use

- The H-Plate 11 Holes is indicated for the internal fixation of fractures, osteotomies, and reconstructive procedures in the craniofacial or mandibular regions. Its H-shaped design offers multidirectional stability and is especially useful in areas requiring broad surface coverage and secure bone fixation.