



CAT. NO.	L-PLATE 3+4 HOLES
706.511L	L-PLATE 3+4 HOLES LEFT THICKNESS 0.6MM
706.511R	L-PLATE 3+4 HOLES RIGHT THICKNESS 0.6MM
707.011L	L-PLATE 3+4 HOLES LEFT THICKNESS 0.7MM
707.011R	L-PLATE 3+4 HOLES RIGHT THICKNESS 0.7MM
707.611L	L-PLATE 3+4 HOLES LEFT THICKNESS 0.8MM
707.611R	L-PLATE 3+4 HOLES RIGHT THICKNESS 0.8MM

Specification

- Plate is available in Titanium.

Indication for use

- The L-Plate (3+4 holes) is used in oral and maxillofacial surgery, particularly for fixation of facial bone fractures, where a right-angle configuration is beneficial. It is commonly used for mandibular angle fractures, orbital rim fractures, or zygomaticomaxillary complex (ZMC) repairs.