




| CAT NO.  | L-PLATE 4+6 HOLES                        |
|----------|------------------------------------------|
| 706.516L | L-PLATE 4+6 HOLES LEFT THICKNESS 0.6MM   |
| 706.516R | L-PLATE 4+6 HOLES RIGHT THICKNESS 0.6 MM |
| 707.016L | L-PLATE 4+6 HOLES LEFT THICKNESS 0.7MM   |
| 707.016R | L-PLATE 4+6 HOLES RIGHT THICKNESS 0.7 MM |
| 707.616L | L-PLATE 4+6 HOLES LEFT THICKNESS 0.8MM   |
| 707.616R | L-PLATE 4+6 HOLES RIGHT THICKNESS 0.8 MM |

### Specification

- Plate is available in Titanium.

### Indication for use

- The L-Plate (4+6 holes) is a larger, L-shaped fixation plate used in maxillofacial surgery to provide rigid stabilization in areas of the facial skeleton where two planes meet at a right angle, such as the mandibular angle or zygomaticomaxillary complex.

