

**LCP MANDIBLE PLATE, ANGLED
BROAD CENTRE, 3 + 3 HOLES**



CAT NO.	LCP MANDIBLE PLATE, ANGLED BROAD CENTRE, 3+3 HOLES
709.356	LCP MANDIBLE Plate, angled broad centre, 3+3 holesTHICKNESS 2 MM

Specification

- Plate is available in Titanium.

Indication for use

- The LCP Mandible Plate (Angled Broad Centre, 3+3 Holes) is designed for stabilizing fractures in the mandible, particularly in the angle and body regions. Its broad central section provides enhanced support, while the angled configuration allows for optimal alignment and fixation, promoting effective healing.

