

## LCP MANDIBLE PLATE, STRAIGHT, WITH CENTRE SPACE, 2+2 HOLES



CAT NO.	LCP MANDIBLE PLATE, STRAIGHT, WITH CENTRE SPACE, 2+2 HOLES
708.511	LCP MANDIBLE Plate, straight, with centre space, 2+2 holes, Thickness 1.25 MM

### Specification

- Plate is available in Titanium.

### Indication for use

- The LCP Mandible Plate with a straight configuration and 2+2 Hole arrangement is designed for the internal fixation of mandibular fractures, particularly in areas requiring stable fixation with minimal bone contact. Its design allows for adaptable contouring to fit the anatomical contours of the mandible, ensuring optimal stabilization during the healing process.