



CAT. NO.	LCP MANDIBLE ADAPTION PLATE
708.064	LCP MANDIBLE Adaption Plate, 4 holes, Thickness 1.0 MM
708.066	LCP MANDIBLE Adaption Plate, 6 holes, Thickness 1.0 MM
708.072	LCP MANDIBLE Adaption Plate, 12 holes, Thickness 1.0 MM
708.080	LCP MANDIBLE Adaption Plate, 20 holes, Thickness 1.0 MM
708.542	LCP MANDIBLE Adaption Plate, 12 holes, Thickness 1.25 MM

## Specification

- Plate is available in Titanium.

## Indication for use

- The LCP Mandible Adaption Plate is utilized in oral and maxillofacial surgery to stabilize fractures of the mandible, particularly in the sub condylar and condylar regions. Its locking mechanism allows for secure fixation without requiring precise anatomical contouring, making it suitable for complex or comminuted fractures

