



CAT . NO.	LCP MANDIBLE DCP PLATE, ANGLED 3+3 HOLES
708.531	LCP MANDIBLE DCP Plate, angled 3+3 holes, Thickness 1.25 MM

Specification

- Plate is available in Titanium.

Indication for use

- The LCP Mandible DCP Plate (Angled 3+3 holes) is designed for the internal fixation of mandibular fractures, particularly in the angle, body, and symphysis regions. Its dynamic compression mechanism allows for controlled axial compression across the fracture site, promoting optimal bone healing.