



CAT NO.	LCP MANDIBLE DCP PLATE
708.506	LCP MANDIBLE DCP Plate, 2+2 holes THICKNESS 1.25MM
708.521	LCP MANDIBLE DCP Plate 3+3 holes THICKNESS 1.25MM
708.872	LCP MANDIBLE DCP Plate, 2+2 holes THICKNESS 1.5 MM
708.876	LCP MANDIBLE DCP Plate, 3+3 holes THICKNESS 1.5 MM

## Specification

- Plate is available in Titanium.

## Indication for use

- The LCP Mandible DCP Plate is used to stabilize fractures in the mandible, especially in the angle, body, and symphysis regions. Its design promotes bone healing by allowing controlled compression and preserving blood supply to the bone.