



CAT NO.	LCP MANDIBLE PLATE, CRESCENT 2+2 HOLES
708.516	LCP MANDIBLE Plate, crescent 2+2 holes, Thickness 1.25 MM

Specification

- Plate is available in Titanium.

Indication for use

- The LCP Mandible Plate (Crescent, 2+2 Holes) is designed for stabilizing fractures in the mandible, particularly in the angle, body, and Para symphysis regions. Its crescent-shaped cross-section provides optimal strength and minimal bone contact, promoting effective healing.