



CAT NO.	LCP MANDIBLE PLATE, CRESCENT, 3+3 HOLES
708.863	LCP MANDIBLE Plate, crescent, 3+3 holes THICKNESS 1.5MM

Specification

- Plate is available in Titanium.

.

Indication for use

- The LCP Mandible Plate (Crescent, 3+3 Holes) is intended for the internal fixation of mandibular fractures, particularly in the angle and body regions. Its crescent-shaped cross-section provides optimal strength and minimal bone contact, promoting effective healing.

