



CAT NO.	LCP MANDIBLE PLATE, STRAIGHT 3+3 HOLES
708.843	LCP MANDIBLE Plate, straight 3+3 holes, THICKNESS 1.5 MM

Specification

- Plate is available in Titanium.

Indication for use

- The LCP Mandible Plate with a straight configuration and 3+3-hole arrangement is designed for the internal fixation of mandibular fractures, particularly in areas requiring stable fixation with minimal bone contact. Its design allows for adaptable contouring to fit the anatomical contours of the mandible, ensuring optimal stabilization during the healing process.