

**LCP MANDIBLE STRUT PLATE,  
8 HOLES, MALLEABLE**



CAT NO.	LCP MANDIBLE STRUT PLATE, 8 HOLES, MALLEABLE
708.098	LCP MANDIBLE Strut Plate, 8 holes, Malleable Thickness 1.0 MM

**Specification**

- Plate is available in Titanium.

**Indication for use**

- The LCP Mandible Strut Plate with 8 holes is a malleable, three-dimensional implant designed to provide enhanced torsional stability in the fixation of mandibular angle fractures. Its curved and buttressed design allows for precise anatomical adaptation, improving resistance to multidirectional forces and promoting optimal healing.