

## LCP MANDIBLE MINI TENSION BAND PLATE, BROAD CENTRE



CAT. NO.	LCP MANDIBLE MINI TENSION BAND PLATE, BROAD CENTRE
708.031	LCP MANDIBLE Mini Tension Band Plate, broad centre, 2+2 holes, malleable
708.036	LCP MANDIBLE Mini Tension Band Plate, broad centre, 2+3 holes, malleable

### Specification

- Plate is available in Titanium.

### Indication for use

- The LCP Mandible Mini Tension Band Plate (Broad Centre) is indicated for fixation of mandibular fractures, especially in the anterior region requiring wider support. It provides stable fixation in symphysis and parasymphysis fractures and is suitable for use in load-sharing and reconstructive procedures.