

# LCP MANDIBLE MINI TENSION BAND PLATE, NARROW CENTRE, 2 + 2 HOLES



<b>CAT NO.</b>	<b>LCP MANDIBLE MINI TENSION BAND PLATE, NARROW CENTRE, HOLES</b>
708.011	LCP MANDIBLE Mini Tension Band Plate, narrow centre, 2+2 holes, 1.0MM THICKNESS PLATE

## Specification

- Plate is available in Titanium.

## Indication for use

- The LCP Mandible Mini Tension Band Plate (Narrow Centre, 2+2 Holes, Malleable) is designed for stabilizing fractures in the mandible, especially in the parasymphysis and body regions. Its malleable structure allows for precise contouring to the bone, while the 2+2 Holes configuration ensures optimal screw placement for stable fixation. The narrow centre design minimizes soft tissue interference, promoting faster healing and reduced risk of complications.