

LCP MANDIBLE MINI TENSION BAND PLATE, NARROW CENTRE



CAT. NO.	LCP MANDIBLE MINI TENSION BAND PLATE,NARROW CENTRE
708.021	LCP MANDIBLE Mini Tension Band Plate, narrow centre, 3+3 holes, malleable Thickness 1.0 MM
708.026	LCP MANDIBLE Mini Tension Band Plate, narrow centre, 3+3 holes Thickness 1.0 MM

Specification

- Plate is available in Titanium.

.

Indication for use

- The LCP Mandible Mini Tension Band Plate (Narrow Centre) is indicated for the fixation of mandibular fractures, particularly in the symphysis and parasymphysis regions. It serves as a tension band for stable fixation in load-bearing and reconstructive procedures.