

LCP MANDIBLE MINI TENSION BAND PLATE BROAD CENTRE, 3 + 3 HOLES



CAT. NO.	LCP MANDIBLE MINI TENSION BAND PLATE BROAD CENTRE, 3+3 HOLES
708.040	LCP MANDIBLE Mini Tension Band Plate broad centre, 3+3 holes, malleable

Specification

- Plate is available in Titanium.

Indication for use

- The LCP Mandible Mini Tension Band Plate (Broad Centre, 3+3 holes) is indicated for the treatment of mandibular fractures in the symphysis and Para symphysis regions. The broad center design offers enhanced support across the fracture site, and the 3+3 hole configuration allows for secure fixation using \varnothing 2.0 mm locking or cortex screws.