

LCP MANDIBLE MINI TENSION BAND PLATE, NARROW CENTRE, 2+2 HOLES, MALLEABLE



CAT. NO.	LCP MANDIBLE MINI TENSION BAND PLATE, NARROW CENTRE, 2+2 HOLES, MALLEABLE
708.006	LCP MANDIBLE Mini Tension Band Plate, narrow centre, 2+2 holes, malleable 1.0MM THICKNESS PLATE

Specification

- Plate is available in Titanium.

Indication for use

• The LCP Mandible Mini Tension Band Plate (Narrow Centre, 2+2 Holes, Malleable) is designed for the internal fixation of mandibular fractures, particularly in the angle, body, and symphysis regions. Its malleable structure allows for precise contouring to the bone, while the 2+2 Holes configuration provides optimal screw placement for stable fixation. The narrow centre design ensures minimal soft tissue interference, promoting faster healing and reduced risk of complications.