

# LCP MANDIBLE TENSION BAND PLATE, PREBENT, 2+2 HOLES



CAT. NO.	LCP MANDIBLE TENSION BAND PLATE, PREBENT, 2+2 HOLES
708.016R	LCP MANDIBLE Tension Band Plate, Prebent, Right, 2+2 holes thickness 1.0 MM
708.016L	LCP MANDIBLE Tension Band Plate, Prebent, Left, 2+2 holes thickness 1.0 MM

## Specification

- Plate is available in Titanium.

## Indication for use

- The LCP Mandible Tension Band Plate (prebent, 2+2 holes, malleable) is designed for the internal fixation of mandibular fractures, particularly in areas subjected to tensile forces. Its prebent configuration and malleability facilitate optimal anatomical adaptation, providing enhanced stability during the healing process.