



CAT. NO.	MIDFACE ADAPTION PLATE 6 HOLES	
706.556	MIDFACE ADAPTION PLATE 6 HOLES	THICKNESS 0.6 MM
707.036	MIDFACE ADAPTION PLATE 6 HOLES	THICKNESS 0.7 MM
707.686	MIDFACE ADAPTION PLATE 6 HOLES	THICKNESS 0.8 MM

Specification

- Plate is available in Titanium.

Indication for use

- The Midface Adaption Plate (6 Holes) is indicated for the stabilization and fixation of fractures, osteotomies, and reconstructive procedures in the midfacial skeleton, including the maxilla and zygomatic regions. It is designed to conform to complex anatomical contours, providing reliable support in delicate facial structures.