



CAT. NO.	MIDFACE ADAPTION PLATE 8 HOLES
706.558	MIDFACE ADAPTION PLATE 8 HOLES THICKNESS 0.6 MM
707.038	MIDFACE ADAPTION PLATE 8 HOLES THICKNESS 0.7 MM
707.688	MIDFACE ADAPTION PLATE 8 HOLES THICKNESS 0.8 MM

Specification

- Plate is available in Titanium.

Indication for use

- The Midface Adaption Plate (8 Holes) is indicated for internal fixation of midfacial fractures, osteotomies, and reconstructive procedures involving areas such as the maxilla and zygomatic bones. Its 8-hole design allows for versatile screw placement and stable adaptation to complex facial anatomy.