



<b>CAT. NO.</b>	<b>Y-PLATE 3 HOLES</b>
706.611	Y-PLATE 3 HOLES THICKNESS 0.6 MM
707.091	Y-PLATE 3 HOLES THICKNESS 0.7 MM

## Specification

- Plate is available in Titanium.

## Indication for use

- The Y-Plate 3 Holes is indicated for the fixation of fractures, osteotomies, or reconstructions in small or angled anatomical areas, such as the midface or mandibular region. Its Y-shaped design offers stable support at bone junctions or branching fracture lines.