



CAT. NO.	METAL	GREENS HEAD METAL FOR THR
138.004S	SS	GREENS HEAD METAL FOR THR 22 /+0 STERILISED
138.005S	SS	GREENS HEAD METAL FOR THR 22 /+7 STERILISED
138.006S	SS	GREENS HEAD METAL FOR THR 22 /+4 STERILISED
138.012S	SS	GREENS HEAD METAL FOR THR 28/+15.5 STERILISED
138.013S	SS	GREENS HEAD METAL FOR THR 28 /+12 STERILISED
138.014S	SS	GREENS HEAD METAL FOR THR 28 /+8.5STERILISED
138.015S	SS	GREENS HEAD METAL FOR THR 28 /+5 STERILISED
138.016S	SS	GREENS HEAD METAL FOR THR 28 /+1.5 STERILISED
138.021S	SS	GREENS HEAD METAL FOR THR 32 /+17 STERILISED
138.022S	SS	GREENS HEAD METAL FOR THR 32 /+13 STERILISED
138.023S	SS	GREENS HEAD METAL FOR THR 32 /+9 STERILISED
138.024S	SS	GREENS HEAD METAL FOR THR 32 /+5 STERILISED
138.025S	SS	GREENS HEAD METAL FOR THR32 /+1 STERILISED
138.031S	SS	GREENS HEAD METAL FOR THR 36 /+12 STERILISED
138.032S	SS	GREENS HEAD METAL FOR THR 36 /+8.5 STERILISED
138.033S	SS	GREENS HEAD METAL FOR THR 36 /+5 STERILISED
138.034S	SS	GREENS HEAD METAL FOR THR 36 /+1.5 STERILISED
138.035S	SS	GREENS HEAD METAL FOR THR 36 /2 STERILISED

## Specification

- GREENS Head metal for THR Available in Stainless Steel.

## Indication for use

- A Metal or ceramic ball is placed on the upper part of the stem. This ball replaces the damaged femoral Head that was removed.

