



CAT. NO.	METAL	UNCEMENTED THR
750.616S	TT	GREENS UNCEMENTED STEM S-6 STERLISED
750.618S	TT	GREENS UNCEMENTED STEM S-8 STERLISED
750.619S	TT	GREENS UNCEMENTED STEM S-9 STERLISED
750.620S	TT	GREENS UNCEMENTED STEM S-10 STERLISED
750.621S	TT	GREENS UNCEMENTED STEM S-11 STERLISED
750.622S	TT	GREENS UNCEMENTED STEM S-12 STERLISED
750.623S	TT	GREENS UNCEMENTED STEM S-13 STERLISED
750.624S	TT	GREENS UNCEMENTED STEM S-14 STERLISED

## Specification

- Greens Uncemented THR available Stainless Steel.

## Indication for use

- The uncemented femoral stem is press fit, double tapered design with a hydroxyapatite (HA) coating and is intended to be implanted in prepared femur medullary canal without bone cement.

