


ANTERIOR CERVICAL PLATE SYSTEM		
CAT NO.	PRODUCT NAME	QTY
	VARIABLE ANGLE CERVICAL PLATE	
690.217	VARIABLE ANGLE CERVICAL PLATE 17.0MM	1
690.219	VARIABLE ANGLE CERVICAL PLATE 19.0MM	1
690.221	VARIABLE ANGLE CERVICAL PLATE 21.0MM	1
690.223	VARIABLE ANGLE CERVICAL PLATE 23.0MM	1
690.225	VARIABLE ANGLE CERVICAL PLATE 25.0MM	1
690.227	VARIABLE ANGLE CERVICAL PLATE 27.5MM	1
690.230	VARIABLE ANGLE CERVICAL PLATE 30.0MM	1
690.232	VARIABLE ANGLE CERVICAL PLATE 32.5MM	1
690.235	VARIABLE ANGLE CERVICAL PLATE 35.0MM	1
690.237	VARIABLE ANGLE CERVICAL PLATE 37.5MM	1
690.240	VARIABLE ANGLE CERVICAL PLATE 40.0MM	1
690.242	VARIABLE ANGLE CERVICAL PLATE 42.5MM	1
690.245	VARIABLE ANGLE CERVICAL PLATE 45.0MM	1
690.247	VARIABLE ANGLE CERVICAL PLATE 47.5MM	1
690.250	VARIABLE ANGLE CERVICAL PLATE 50.0MM	1
690.255	VARIABLE ANGLE CERVICAL PLATE 55.0MM	1
690.260	VARIABLE ANGLE CERVICAL PLATE 60.0MM	1
	LCP VARIABLE ANGLE CERVICAL SCREW 3.5MM	
690.312	LCP VARIABLE ANGLE CERVICAL SCREW 3.5MM X 12.0MM	8
690.314	LCP VARIABLE ANGLE CERVICAL SCREW 3.5MM X 14.0MM	8
690.316	LCP VARIABLE ANGLE CERVICAL SCREW 3.5MM X 16.0MM	8
690.318	LCP VARIABLE ANGLE CERVICAL SCREW 3.5MM X 18.0MM	8
690.320	LCP VARIABLE ANGLE CERVICAL SCREW 3.5MM X 20.0MM	8
	LCP VARIABLE ANGLE CERVICAL RESCUE SCREW 3.8 FOR CERVICAL PLATE DESIGN TYPE-1 (TITANIUM)	
690.332	LCP VARIABLE ANGLE CERVICAL RESCUE SCREW 3.8 / 12MM FOR CERVICAL PLATE DESIGN TYPE-1 (TITANIUM)	8
690.333	LCP VARIABLE ANGLE CERVICAL RESCUE SCREW 3.8 / 14MM FOR CERVICAL PLATE DESIGN TYPE-1 (TITANIUM)	8
690.334	LCP VARIABLE ANGLE CERVICAL RESCUE SCREW 3.8 / 16MM FOR CERVICAL PLATE DESIGN TYPE-1 (TITANIUM)	8



ANTERIOR CERVICAL PLATE INSTRUMENT SET		
SR NO	NAME OF INSTRUMENTS	QTY
1	PLATE HOLDER-I FOR CERVICAL PLATE.	1
2	SCREW HOLDER FOR CERVICAL SCREW	1
3	LOCKING PIN SCREW DRIVER	1
4	DISTRACTOR PIN TIGHTENER	1
5	HEX SCREW DRIVER SW2.5	1
6	MULTIAXIAL SCREW DRIVER SW2.5	1
7	QUICK COUPLING HANDLE FOR CERVICAL PLATE	1
8	DISTRACTOR PIN FOR CERVICAL	1
9	IMPLANTER OF LOCATION SCREW FOR CERVICAL PLATE	1
10	LIMITATOR WRENCH FOR CERVICAL PLATE	1
11	QUICK COUPLING HEX SCREW DRIVER SW2.5 FOR CERVICAL PLATE	1
12	DRILL BIT 2.2MM FOR CERVICAL LATE	1
13	LIMITED DRILL BIT 2.2MM FOR CERVICAL PLATE	1
14	IMPLANTS CASE FOR CERVICAL PLATE	1
15	DEPTH GOUGE FOR CERVICAL PLATE	1
16	CERVICAL LOACK SPANNER	1
17	BONE AWL FOR CERVICAL	1
18	MULTIAXIAL DRILL GUIDER (B2.2MM	1
19	DRILL GUIDER (B2.2MM	1
20	PLATE HOLDER - 2 FOR CERVICAL PLATE	1
21	PLATE BENDER WRENCH	1
22	CERVICAL DISTRACTOR	1



Specification

- All implants are available in TT.
- All instruments are available in SS.

Indication for use

- The Anterior Cervical Instrument Set is indicated for use in anterior cervical spine surgeries, such as discectomy, fusion, and decompression procedures. It is used to prepare the intervertebral space, insert implants, and provide stabilization for cervical spine reconstruction.