



BROKEN SCREW REMOVAL INSTRUMENT SET			
SR. NO.	NAME OF INSTRUMENT	SIZE (MM)	QTY
1	Square Screwdriver I	φ1.2	1
2	Square Screwdriver II	φ1.5	1
3	T-Shaped Wrench		1
4	Ring drill I	Large	1
5	Ring drill II	Middle	1
6	Ring drill III	Small	1
7	Cross screwdriver		1
8	Flat screwdriver I	Small	1
9	Flat screwdriver II	Large	1
10	Torx screwdriver I	T8	1
11	Torx screwdriver II	T9	1
12	Torx screwdriver III	T10	1
13	Torx screwdriver IV	T15	1
14	Torx screwdriver V	T18	1
15	Torx screwdriver VI	T20	1
16	Torx screwdriver VII	T25	1
17	Torx screwdriver VIII	T30	1
18	Hex screwdriver I	SW1.5	1
19	Hex screwdriver II	SW2.0	1
20	Hex screwdriver III	SW2.5	1

21	Hex screwdriver IV	SW3.0	1
22	Hex screwdriver V	SW3.5	1
23	Hex screwdriver VI	SW4.0	1
24	Hex screwdriver VII	SW4.5	1
25	Hex screwdriver VIII	SW5.0	1
26	Hex screwdriver IX	SW5.5	1
27	Hex extractor I	1.2mm	1
28	Hex extractor II	2.5mm	1
29	Hex extractor III	2.7mm	1
30	Hex extractor IV	3.5mm	1
31	Hex extractor V	4.0mm	1
32	Broken screw extractor I	φ2.0	1
33	Broken screw extractor II	φ2.7	1
34	Broken screw extractor III	φ4.0	1
35	Broken screw extractor IV	φ5.0	1
36	Broken screw extractor V	φ6.5	1
37	Broken screw extractor VI	φ8.0	1
38	Reamer		1
39	Sharp hook		1
40	T-shaped quick change handle		1
41	Straight quick change handle		1
42	Broken screw Puller I	φ4.5	1
43	Broken screw Puller II	φ3.5	1
44	Screw holder		1
45	Round bone knife	8mm	1
46	Drill I	φ5	1
47	Drill II	φ6	1
48	Drill III	φ7	1

Specification

- All instruments are available in SS.

Indication for use

- The Broken Screw Removal Instrument Set is indicated for extraction of fractured or stripped screws from bone during orthopaedic or trauma surgeries. It facilitates safe and efficient removal without damaging the surrounding bone or implant site.

