




ELASTIC NAIL (TEN S) SYSTEM	
CAT.NO.	PRODUCTNAME
	ELASTIC NAIL (TEN S)
673.015	TITANIUM ELASTIC NAIL (TEN) - 1.5 MM
673.020	TITANIUM ELASTIC NAIL (TEN) - 2.0 MM
673.025	TITANIUM ELASTIC NAIL (TEN) - 2.5 MM
673.030	TITANIUM ELASTIC NAIL (TEN)- 3.0 MM
673.035	TITANIUM ELASTIC NAIL (TEN) - 3.5 MM
673.040	TITANIUM ELASTIC NAIL (TEN) - 4.0 MM
673.045	TITANIUM ELASTIC NAIL (TEN) - 4.5 MM

ELASTIC NAIL (TENS) INSTRUMENT SET

SR. NO.	NAME OF INSTRUMENT	QTY
1	TENSE NAIL RAM ROD	1
2	TENSE NAILAWL	1
3	TENSE NAIL IMPACTOR HAMMER	1
4	TENSE NAIL EXTRACTOR HAMMER	1
5	TENSE NAIL PROFILE BENDER	3+1
6	TENSE NAIL EXTRACTOR WITH CHUCK	1
7	TENSE NAIL IMPACTOR 1.5 MM X 5.0 STRAIGHT	1
8	TENSE NAIL IMPACTOR 1.5 MM X 2.0 STRAIGHT	1
9	TENSE NAIL IMPACTOR 1.5 MM X 5.0 BEVELLED	1

Specification

- Implants are available in TT.
- Instruments are available in SS.

Indication for use

• The Titanium Elastic Nail System is used for the treatment of diaphyseal and metaphyseal fractures in long bones. It is primarily indicated for paediatric fractures of the femur, tibia, humerus, and forearm bones. The system provides stable, flexible fixation allowing early mobilization and growth plate preservation. It is ideal for minimally invasive internal fixation in skeletally immature patients.