



I.L. EXPERT TIBIA NAIL SYSTEM

I.L. EXPERT TIBIAL NAIL Ø 8 MM

I.L. EXPERT TIBIAL NAIL Ø 9 MM

I.L. EXPERT TIBIAL NAIL Ø 10 MM

I.L. EXPERT TIBIAL NAIL Ø 11 MM

LOCKING BOLT (SELF TAPPING) Ø 4.9 MM

LOCKING BOLT (SELF TAPPING) Ø 3.9 MM

END CAP FOR NAIL EXTENSION

I.L. EXPERT TIBIA NAIL INSTRUMENT SET		
CAT.NO.	NAME OF INSTRUMENTS	QTY
903.501	EXPERT TIBIA NAIL JIG	1
903.105	INSERTION DRIVING HEAD	1
903.110	RAM	1
903.115	RAM ROD	1
903.502	NAIL CONNECTING BOLT	1
903.503	PROTECTION SLEEVE 10.0 X 8.0 MM	1
903.504	DRILL SLEEVE 8.0 X 4.0 MM	1
903.505	TROCAR 8.0 MM	1
903.120	WRENCH	1
903.145	TISSUE PROTECTOR	1
903.150	DEPTH GAUGE FOR INTERLOCKING 4.9 MM	1
903.160	GUIDE ROD 2.5 MM X 95 CM	1
903.155	GUIDE WIRE HOLDER	1
439.500	SCREW DRIVER (HEXAGONAL) 4.5 MM FOR I.L.	1
903.165	REAMING ROD WITH OLIVE POINT 2.5 MM X 95 CM	1
426.140	DRILL BIT 4.0 MM X 225 MM NON-COUPLING	1
401.100	KUNTSCHER'S DIAMOND POINTED AWL CURVED	1
903.205	WIRE SLEEVE 8 MM X 2.0 MM	1
903.225	TEFLON TUBE	1
903.105	DRIVING HEAD	1

Specification

- Implants are available in TT & SS.
- Instruments are available in SS.

Indication for use

- The Expert Tibia Nail System is indicated for intramedullary nailing of tibial fractures, including proximal, shaft, and distal tibial fractures. It is used to prepare the medullary canal, insert the tibial nail, and secure fixation for stable fracture management and early mobilization.